

# MultiView™ Cleaning and Care Guide

MPI utilizes high quality materials on all products, including 4-Way Stretch Knit by Sioen Performance NV.

## 4-way stretch knit, cast transfer coated urethane

Internal testing of commonly utilized products resulted in the following recommendations:

### Disinfectant Use:



The information published in this care guide refers to the performance of Sioen Performance NV products in specific tests conducted under laboratory conditions. Results may vary under actual conditions.

- SANI Professional disinfectant wipes
- Formula 409® All-Purpose spray cleaner
- Isopropyl Alcohol
- Bleach
- Quaternary Ammonium Solution
- Virucide
- Water

## Cleaning After Disinfection:



The following cleaners should be used with soft cloth or damp sponge to remove chemical residue.

- Formula 409® All-Purpose spray cleaner
- Mild liquid soap and water

### Remove most stains with one of the following steps:

**Step 1:** eye shadow, mascara, grease, black felt tip pen, crayon, chocolate, suntan lotion and hair/body oils\*. Remove excess spill with a damp cloth. Clean with a 1:1 mix of Mild soap and water. Next, rinse with clean water and dry.

**Step 2:** red lipstick  
Use a straight application of concentrated cleaners such as Formula 409® spray cleaner. Wipe with a clean cloth, rinse with water and dry.

**Step 3:** blood, urine, feces, oil base paint, tar and asphalt.  
Use a 1:10 mix of bleach and water. Rinse with clean water and dry.

**Step 4:** All other tough stains  
Use a 1:10 mix of isopropyl alcohol and water. Rinse thoroughly with clean water and pat dry.

**Note:** For cleaning that requires step 4, use a soft cotton cloth saturated with the cleaning material and rub the stain in circles 10 times. Pat dry with another soft cotton cloth and check the results.

**After cleaning, rinse area with clean water, then dry with clean, lint-free cloth. Failure to remove chemical residue will cause discoloration of the coated fabric surface to manifest earlier in the life cycle of the product.**