

# UltraScan Premier Table™

The ultimate multi-modality ultrasound platform

MEDICAL POSITIONING, INC.  
www.medicalpositioning.com  
800-593-3246



ergonomic \* innovative \* proven

Optional features may be shown. Actual color may vary depending on monitor resolution and print quality.

## STANDARD FEATURES

<b>Height Adjustability:</b>	<b>Maximum Patient Capacity:</b>	<b>Fowler (Head Up):</b>	<b>Trendelenburg/Reverse Trendelenburg:</b>	<b>Control System:</b>	<b>Braking:</b>
24" to 34"	500 lbs	90°	15° /25°	Hand Controller	Single Pedal Braking
<b>Stirrups:</b>	<b>Foot Down Section:</b>	<b>Ergonomical Features:</b>	<b>Additional Features:</b>		
Adjustable Hideaway	Manually adjustable / Removable  Built-In Folding Foot Platform	Patient Arm Boards and Back Support  Sonographer Arm Rest  Patient Head Support	Hideaway Safety Rails		

## PRODUCT SPECIFICATIONS

<b>Length:</b>	<b>Width:</b>	<b>Weight:</b>
72"	28"	290 lbs
<b>Foam:</b>	<b>Electrical:</b>	<b>Vinyl:</b>
California Technical Bulletin 117 - Section E	120 VAC, 50/60 Hz, 2.0A Max  230 VAC, 50/60 Hz, 1.0A Max	FMVSS 302  Boston Fire Code - BFD IX-I  NY/NJ Port Authority

## PRODUCT BENEFITS

Perform multiple ultrasound procedures using one table  
 Optimizes patient access and positioning to reduce image acquisition time  
 Maximum lower extremity vessel dilation  
 Protects sonographers from potential career-ending injuries by improving ergonomics  
 Bariatric capable

## CERTIFICATIONS

<b>Industry Standards:</b>	<b>Warranty:</b>
UL 60601-1, CAN/CSA C22.2 No. 601.1	1 Year

## OPTIONAL FEATURES

Foot Controller  
 IV Pole and Holder