CathErgometer

PRODUCT FEATURES (*SOME OPTIONAL FEATURES SHOWN)

Cath Ergometer[™]







Cath **Ergometer**

Standard Features

Removable Ergometer to Specific Cath or Nuclear Table of choice

- Patient Coach with Audible and Visual Patient Guidance
- Ergometer Boots
- Ergometer Controller
 - 10 Pre-Programmed, 4 Programmable, 1 Manual Protocol

Optional Features

Storage Cart

Ergometer Interface

Additional Product Information

Protocols:

10 Pre-Programmed, 4 Programmable, 1 Manual

Electrical:

120 VAC, 50/60 Hz Max Watts 300

Warranty:

2 Years

Industry Standards:

UL 60601-1, CAN/CSA C22.2 No. 601.1)

