


Stress EchoBed® Ergometer Replacement

medicalpositioning.com
800-593-3246

PRODUCT FEATURES (*SOME OPTIONAL FEATURES SHOWN)

Stress EchoBed® Ergometer Replacement



Stress EchoBed® Ergometer Replacement

STANDARD FEATURES

Adjustable/Removable Ergometer

- Patient Coach with Audible and Visual Patient Guidance
- Ergometer Boots
- Ergometer Controller
 - 10 Pre-Programmed, 4 Programmable, 1 Manual Protocol

OPTIONAL FEATURES

Storage Cart

Ergometer Interface

ADDITIONAL PRODUCT INFORMATION

Electrical:

120 VAC, 50/60 Hz

Max Watts 300

Protocols:

10 Pre-Programmed, 4 Programmable, 1 Manual

Warranty:

2 Years

Industry Standards:

UL 60601-1, CAN/CSA C22.2 No. 601.1

