

Stress EchoBed®



800.593.3246 • info@medicalpositioning.com

STANDARD FEATURES

Height Range: 24"-34"	Load Capacity: 500 lbs.	Patient Lateral Tilt Left: 40°
Trendelenburg: +15°/-15° Trendelenburg Positioning	Ergonomic Features: Dual Rapid Release Imaging Drop Sections	Braking: Individual Locking Casters

ADDITIONAL FEATURES



Removable Safety Handrails



Positioning Wedge



Patient Safety Belt



Patient Shoulder Support



Removable Ergometer

OPTIONAL FEATURES



Ergometer Storage Cart



Ergometer Track Cover




Sonographer Extension



Articulating Padded Arm Board



IV Pole and Holder



Pediatric/Geriatic Adapter



Paper Roll Holder and Cutter



Grab Handle

ADDITIONAL PRODUCT INFORMATION

Table Length (without Ergometer) 72"	Table Length (with Ergometer) 81"-94"
Product Width: 30"	Product Weight: 492 lbs. with Ergometer 420 lbs. without Ergometer
Control System: Hand Controller	Warranty: 5 Years
Standard Color: Pearl	




Stress EchoBed®



MedicalPositioning

SUPINE BIKE STRESS
ECHOCARDIOGRAPHY



OPTIMAL PATIENT POSITIONING

A powerful solution for stress echocardiograms, the Stress EchoBed® enables clinicians to capture images during peak exercise—even in patients who cannot use a treadmill.

Stress EchoBed[®]

BENEFITS AND FEATURES

- Induces symptoms at lower heart rates
- Provides imaging during intermediate stages to identify the ischemic threshold
- Improves diagnostic capabilities for valvular pathology and pulmonary vascular problems
- Eliminates rapidly falling heart rates seen post-treadmill
- Avoids the need for post-exercise transfer
- Evaluates diastolic pressures and manages heart failure patients
- 4 point safety system
- Bariatric capable

