

— UltraMamm™

MEDICAL POSITIONING, INC.

www.medicalpositioning.com

800-593-3246



Convertible platform for 360° access to breast

Chair to table positioning for multi use

The UltraMamm maximizes department space by incorporating upright and decubitus stereotactic breast biopsy, as well as mammography into one bariatric capable platform, increasing throughput and improving continuity of care.

- Easy transition from chair to table provides an all-encompassing platform for varying patient population
- Accommodates both upright and decubitus positions for stereotactic breast biopsy, as well as mammography and breast ultrasound
- Bariatric capable patient weight capacity of 500 lbs establishes continuity of care
- Reduce turnover time and increase throughput, even in small exam rooms
- Accessories Include: removable arm rests, back support, adjustable and removable head support


Optional features may be shown. Actual color may vary depending on monitor resolution and print quality.

STANDARD FEATURES

Height Adjustability:	Maximum Patient Capacity:	Trendelenburg:	Fowler:	Drop Sections:	Supported Patient Positions:
22" to 38"	500 lbs	15° to -15°	0° to 90°	1 on each side: for clinician and receptor	Upright and Decubitus
Foot Platform:	Control System:	Braking:	Additional Features:		
Electrically adjustable calf section with fold-away foot board	Hand Controller with memory and lock	Single pedal braking	Adjustable and Removable Head Support Adjustable and Removable Back Support Adjustable and Removable Arm Rests Patient arm support sling		

PRODUCT SPECIFICATIONS

Length:	Surface Width:	Length of Wheel Base:	Foam:	Electrical:	Vinyl:
75" with Head Support 64.5" without Head Support	33.5" with arm rests 28.5" without arm rests	33.5"	TB 117-2013 Section E	100-240 VAC, 50/60 Hz 4.5A max	T.B. 117-2013 NFPA 260- Cover Fabric – Class I UFAC Fabric – Class I

Certifications:	Warranty:	Standard Color:	OPTIONAL FEATURES	
UL ANSI/AAMI ES60601-1 (2005) + AMD (2012), CAN/CSA C22.2 No. 60601 (2014)	2 Years		IV Pole and Holder Battery Operation Back Support Wedge Cushion	Mounted Height Foot Control Underbed LED Lighting Instrument Tray