

— UltraMamm™

The ultimate platform for Mammography and Breast Biopsy

MEDICAL POSITIONING, INC.
www.medicalpositioning.com
800-593-3246

Weight capacity up to 500 lbs.

90° Fowler Positioning

Articulating patient arm boards
and patient back supports

Optional $\pm 15^\circ$ Trendelenburg Positioning

Battery operated

22" to 38" height

Electric built-in calf section
and folding foot platform




STANDARD FEATURES

Height Adjustability:	Maximum Patient Capacity:	Fowler (Head Up):	Calf Section:	Braking:
22" to 38"	500 lbs	0 - 90°	0 - 83°	Individual Locking Casters
Control System:	Foot Down Section:	Arm Boards:		
Hand Controller Included Battery Functionality	Electrically Adjustable Built-In Folding Foot Platform	Dual Articulating Patient Arm Boards and Back Support		

PRODUCT SPECIFICATIONS

Length:	Width:	Weight:
72.5"	28.5"	350 lbs
Foam:	Electrical:	Vinyl:
TB 117-2013 Section E	100 - 240 VAC, 50/60 Hz, 4.5A Max	TB 117-2013 NFPA 260 Cover Fabric - Class 1 UFAC Fabric - Class 1

CERTIFICATIONS

Industry Standards:	Warranty:	Standard Color:
UL ANSI/AAMI ES60601-1 (2005) + AMD (2012), CAN/CSA C22.2 No. 60601 (2014)	2 Years	

PRODUCT BENEFITS

- Access lesions in upper/outer quadrant
- For use with any stereotactic add-on system
- Enhance ergonomic posture as technologists work upright
- Bariatric capable

OPTIONAL FEATURES

- Articulating Arm Board Drop Section
- Carotid Head Rest
- Foot Control
- Handrails with Hideaway Brackets
- Head Rail
- IV Pole and Holder
- Paper Roll Holder and Cutter
- Push Handles
- Safety Straps
- Single Pedal Braking +/-15°
- Trendelenburg Positioning
- Underbed LED Lighting