

VasScan™ Table X

medicalpositioning.com
800-593-3246

ADA COMPLIANT 19" LOW HEIGHT TO 38" MAXIMUM HEIGHT

The VasScan™ X is designed to help minimize the risk of repetitive stress injury, reduce exam time and improve results while increasing bariatric capability and meeting the ADA height transfer requirement.

- 19"-38" HEIGHT RANGE MEETS THE AMERICANS WITH DISABILITIES ACT ADA LOW HEIGHT REQUIREMENT OF 17"-19" FOR PATIENT TRANSFER, WHILE ENSURING ERGONOMIC STANDING POSTURE FOR STAFF
- 30° REVERSE TRENDLENBURG IS DESIGNED TO INCREASE VENOUS POOLING AND PROMPT RAPID RECOGNITION AND ACCURATE RESULTS
- THE WEIGHT-BEARING ADJUSTABLE AND REMOVABLE FOOT PLATFORM ACCOMMODATES PATIENTS OF VARYING HEIGHTS WHILE ENSURING PATIENT SAFETY DURING LOWER VENOUS STUDIES
- BARIATRIC CAPABLE PATIENT WEIGHT CAPACITY OF 600 LBS ESTABLISHES CONTINUITY OF CARE
- 65° FOWLER POSITIONING INCREASES COMFORT OF/FOR PATIENTS WITH SHORTNESS OF BREATH
- OPTIONAL IMAGING AND SONOGRAPHER DROP SECTIONS PROVIDE UNHINDERED ACCESS AND REDUCE IMAGE ACQUISITION TIME FOR ECHOCARDIOGRAMS

(SOME OPTIONAL FEATURES SHOWN)



CARDIAC IMAGING
VASCULAR IMAGING



MedicalPositioning



STANDARD FEATURES

Height Range:

19"-38"

Patient Weight Capacity:

600 lbs.

Braking:

Individual Locking Casters

Control System:

Hand Control

Foot Platform:

Manually Adjustable and Removable

Fowler Positioning:

65°

Trendelenburg Positioning:

+15°/-30°

ADDITIONAL FEATURES

- Hand Control Storage
- Foley Bag/Accessory Hook

NEW AND IMPROVED FEATURES

- ADA Compliant transfer height of 19"
- Maximum height of 38" ergonomically accommodates standing staff
- Increased Patient weight capacity to 600 lbs
- Increased warranty to 2 years
- Memory Positioning - Store up to 2 positions with the touch of a button
- Underbed Lights
- Battery operation - avoid cords when preferred
- Barrier free edges with new hideaway brackets
- Updated and durable hand control
- Ability to utilize both hand and foot controls
- Carotid Head Support to ergonomically accommodate all scanning needs

ADDITIONAL PRODUCT INFORMATION

Table Length:

80.5"

Table Width:

28"

Warranty:

2 Years

Electrical:

100-240 VAC,
50/60 Hz, 5A

Vinyl:

Modena EcoSense
from Spradling®

Standard Color:



Industry Standards:

AAMI ES60601-1 (2005),
CAN/CSA C22.2 No. 60601 (2014)
IEC 60601-1: 2005 + CORR. 1:2006 + CORR.
2:2007 + AM1:2012
(or IEC 60601-1: 2012 reprint)

OPTIONAL FEATURES

- Echocardiogram Drop Sections
- Sonographer Extension
- 2-Way Sonographer Drop Section/
Patient Back Support
- Non-Pinch Flaps
- Single Pedal Braking
- Foot Controller (nonmemory)
- Battery Operation
- Safety Handrails w/Hideaway
Brackets
- Fold-Down Safety Handrails
- Safe-T-Wedge
- Pediatric/Geriatric Adapter
- IV Pole and Holder
- Padded Arm Board
- Paper Roll Holder and Cutter
- Carotid Head Support
- Underbed Lights

