


UltraScan Plus™

The ultimate multi-modality ultrasound platform

MEDICAL POSITIONING, INC.
www.medicalpositioning.com
800-593-3246



90° Fowler positioning

2-way sonographer drop section/
Patient back support (optional)

Imaging drop section (optional)

Weight capacity up to 500 lbs

25° Trendelenburg and 15° Reverse Trendelenburg
with built-in folding foot platform, 15°/15° without
(optional)

24" to 34" height

Adjustable hideaway stir ups

Manually adjustable and
removable foot down section



ergonomic * innovative * proven

Optional features may be shown. Actual color may vary depending on monitor resolution and print quality.

STANDARD FEATURES

Height Adjustability:	Maximum Patient Capacity:	Fowler (Head Up):	Foot Down Section:	Stirrups:	Control System:	Braking:
24" to 34"	500 lbs	90°	Manually adjustable / Removable	Adjustable Hideaway	Hand Controller	Individual Locking Casters

PRODUCT SPECIFICATIONS

Length:	Width:	Weight:
72"	28"	290 lbs
Foam:	Electrical:	Vinyl:
California Technical Bulletin 117 - Section E	120 VAC, 50/60 Hz, 2.0A Max 230 VAC, 50/60 Hz, 1.0A Max	FMVSS 302 Boston Fire Code - BFD IX-I NY/NJ Port Authority

CERTIFICATIONS

Industry Standards:	Warranty:
UL 60601-1, CAN/CSA C22.2 No. 601.1	2 Years

PRODUCT BENEFITS

- Optimizes patient access and positioning to reduce image acquisition time
- Maximum lower extremity vessel dilation
- Protects sonographers from potential career-ending injuries by improving ergonomics
- Bariatric capable

OPTIONAL FEATURES

- Drop Sections
 - Imaging Drop Section
 - 2-Way Sonographer Drop Section/Patient Back Support
- Trendelenburg/Reverse Trendelenburg: 25°/15° with Built-In Folding Foot Platform
- Trendelenburg/Reverse Trendelenburg: 15°/15° without Built-In Folding Foot Platform
- Single Pedal Braking
- Safety Rails
- Foot Controller
- Sonographer Extension
- Pediatric/Geriatric Adapter
- IV Pole and Holder
- Padded Arm Board
- Paper Roll Holder and Cutter