




The ultimate platform for video swallow studies

MEDICAL POSITIONING, INC.
www.medicalpositioning.com
800-593-3246



ergonomic * innovative * proven

Optional features may be shown. Actual color may vary depending on monitor resolution and print quality.

STANDARD FEATURES

Height Adjustability:	Maximum Patient Capacity:	Fowler:	Calf Section:	Arm Rests:	Radiolucent Imaging:	Seat Rotation:
19" to 33"	600 lbs	90°	0° to 75°	Fold Away and Height Adjustable	32" x 13"	90° Left and Right
Powered Operation:	Control System:	Braking:	Additional Features:			
Battery Operated	Hand Controller	Single Pedal Braking	Oxygen Tank Holder Patient Safety Straps Transport Handles			

PRODUCT SPECIFICATIONS

Table Length:	Chair Length:	Width:	Weight:
79.25"	33"	26"	365 lbs
Foam:	Electrical:	Vinyl:	
California Technical Bulletin 117 - Section E	100-230 VAC. 50/60 Hz, 5.0A Max	FMVSS 302 Boston Fire Code - BFD IX-I NY/NJ Port Authority	

CERTIFICATIONS

Industry Standards:	Warranty:
UL ANSI/AAMI ES60601-1 (2005) + AMD (2012), CAN/CSA C22.2 No. 60601 (2014)	1 Year

OPTIONAL FEATURES

IV Pole and Holder
Charger and Second Battery