

ADA compliant 19" low height to 38" maximum height

- 19"-38" height range meets the Americans with Disabilities Act (ADA) low height requirement of 17"-19" for patient transfer, while ensuring ergonomic standing posture for staff
- Optional battery operation
- Dual rapid release drop sections provide unhindered access to the patient
- 65° Fowler positioning increases comfort of/for patients with shortness of breath
- +/- 15° Trendelenburg helps restore patient's blood pressure
- Improved features: new handrail design, improved single pedal braking system, Foley bag hook, hand control storage, battery operation



Deepest imaging drop section on the market accommodates more patients, assuring continuity of care in your facility



Sonographer drop section/
patient back support and rapid
release imaging drop section

Sonographer extension right
handed sonographer support
(optional)

65° Fowler positioning
(optional)

Included: Patient
Positioning wedge

15° reverse Trendelenburg
and 15° Trendelenburg
(optional)

Bariatric capable patient
weight capacity up to 600 lbs.

Static load capacity up
to 1,000 lbs.

30" width

Battery (optional)

Single pedal braking

19" to 38" height
at top of mattress

STANDARD FEATURES

Height Adjustability:	Maximum Patient Capacity:	Drop Sections:	Control System:	Braking:
19" to 38"	600 lbs	Imaging Drop Section	Hand Controller w/ Memory Positioning	Single Pedal Braking
	Maximum Load Capacity: 1,000 lbs	2-Way Sonographer Drop Section/Patient Back Support		

Additional Features:

Safety Rails w/
hideaway brackets
Paper Roll Holder and Cutter
Safe-T-Wedge

Hand control storage
Non-Pinch Flaps
Foley Bag/Accessory Hook

PRODUCT SPECIFICATIONS

Length:	Width:
80.5"	30"

Foam:	Electrical:	Vinyl:
California Technical Bulletin 117 - Section E	100-240VAC 50/60Hz 5A	FMVSS 302 Boston Fire Code - BFD IX-I NY/NJ Port Authority

CERTIFICATIONS

Industry Standards:	Warranty:
UL ANSI/AAMI ES60601-1 (2005) + AMD (2012), CAN/CSA C22.2 No. 60601 (2014)	5 Years

NEW AND IMPROVED FEATURES

ADA Compliant transfer
height of 19"

Maximum height of 38"
ergonomically accommodates
standing sonographers

Increased Patient weight
capacity to 600 lbs

Battery operation -
avoid cords when preferred

Barrier free edges with
new hideaway brackets

Updated and durable hand control

Ability to utilize both hand
and foot controls

Carotid Head Support to
ergonomically accommodate
all scanning needs

Underbed Lights

OPTIONAL FEATURES

Sonographer Extension

Trendelenburg/Reverse Trendelenburg: 15°/15°
-stabilizes patient blood pressure

Fowler (Head Up) Positioning: 65°
-increases comfort of patients with shortness of breath

Battery Operation up to a full day of operation with
no need for plugging in

Foot Controller (non-memory)

Pediatric/Geriatric Adapter

IV Pole and Holder

Padded Arm Board

Carotid Head Support

Underbed Lights