

UltraScan Versa™ Premier

MEDICAL POSITIONING, INC.
www.medicalpositioning.com
800-593-3246

The ultimate multi-modality ultrasound platform

Weight capacity up to 500 lbs

90° Fowler positioning

Patient back support

Articulating patient arm boards

15° Trendelenburg and 25° Reverse Trendelenburg

Adjustable hideaway stirrups

22" to 38" height

Electric built-in calf section and folding foot platform

Single pedal braking



Optional features may be shown. Actual color may vary depending on monitor resolution and print quality.

UltraScan Versa™ Premier Table

STANDARD FEATURES

Height Adjustability:	Maximum Patient Capacity:	Fowler (Head Up):	Trendelenburg/Reverse Trendelenburg:	Calf Section:	Braking:
22" to 38"	500 lbs	0 - 90°	15° /25°	0 - 83°	Single Pedal Braking
Control System:	Foot Down Section:	Stirrups:	Handrails:	Arm Boards:	
Hand Controller	Electrically adjustable Built-In Folding Foot Platform	Adjustable and Hideaway	Removable with Hideaway Bracket	Dual Articulating Patient Arm Boards and Back Support	


PRODUCT SPECIFICATIONS

Length:	Width:	Weight:
72.5"	28.5"	350 lbs
Foam:	Electrical:	Vinyl:
TB 117-2013 Section E	100 - 240 VAC, 50/60 Hz, 4.5A Max	TB 117-2013 NFPA 260 Cover Fabric - Class 1 UFAC Fabric - Class 1

PRODUCT BENEFITS

- Perform multiple ultrasound procedures using one table
- Optimize patient access and positioning to help reduce image acquisition time
- Help protect sonographers from potential career-ending injuries by improving ergonomics
- Bariatric capable
- Increase lower extremity vessel dilation

CERTIFICATIONS

Industry Standards:	Warranty:	Standard Color:
UL ANSI/AAMI ES60601-1 (2005) + AMD (2012), CAN/CSA C22.2 No. 60601 (2014)	5 Years	

OPTIONAL FEATURES

- Carotid Head Rest
- Foot Control
- Head Rail
- IV Pole and Holder
- Paper Roll Holder and Cutter
- Push Handles
- Safety Straps
- Underbed LED Lighting