

— DBI Table™

The ultimate platform for Decubitus Stereotactic Breast Biopsy

MEDICAL POSITIONING, INC.
www.medicalpositioning.com
800-593-3246



Interchangeable head and foot sections

Adjustable/removable Back support

Interchangeable head and foot sections

Weight capacity up to 500 lbs

24" to 36" height

Four imaging drop sections to accommodate more patients without the need for room turnover

Two drop sections on each side to allow space for the image receptor and physician



ergonomic * innovative * proven

Optional features may be shown. Actual color may vary depending on monitor resolution and print quality.

STANDARD FEATURES

Height Adjustability:	Maximum Patient Capacity:	Drop Sections:	Control System:	Braking:
24" to 36"	500 lbs	Four imaging drop sections; two on each side to allow space for the image receptor and physician	Hand Controller	Individual Locking Casters

Additional Features:

Back Support
Interchangeable Head and Foot Sections

PRODUCT SPECIFICATIONS

Length:	Width:	Weight:
75"	27"	315 lbs
Foam:	Electrical:	Vinyl:
California Technical Bulletin 117 - Section E	120 VAC, 50/60 Hz, 1.6A Max 230 VAC, 50/60 Hz, 1.0A Max	California T.B. 117-2013 NFPA 260- Cover Fabric - Class I UFAC Fabric - Class I

CERTIFICATIONS

Industry Standards:	Warranty:	Standard Color:
UL 60601-1, CAN/CSA C22.2 No. 601.1	1 Year	

PRODUCT BENEFITS

- Increase ergonomic access to breast in decubitus positioning, including lesions in upper/outer quadrant
- Compatible with any stereotactic add-on system
- Drop sections provide close access for image receptor
- Enhance ergonomic posture as technicians work upright in drop sections
- Back support allows the patient to take deep breaths while keeping the chest wall in contact with the image receptor
- Minimizes risk of falling associated with fainting during an upright procedure
- Bariatric capable

OPTIONAL FEATURES

Foot Controller