

EchoTable™

The ultimate platform for Cardiac Ultrasound

MEDICAL POSITIONING, INC.
www.medicalpositioning.com
800-593-3246



STANDARD FEATURES

Height Adjustability:	Maximum Patient Capacity:	Drop Sections:	Control System:	Braking:	Ergonomic Features:
24" to 34"	500 lbs	Imaging Drop Section Sonographer Drop Section	Hand Controller	Individual Locking Casters	2-Way Sonographer Drop Section/ Patient Back Support
Additional Features:					

Paper Roll Holder and Cutter

Imaging Drop Section and Back Support

PRODUCT SPECIFICATIONS

Length:	Width:	Weight:
79.5"	28"	300 lbs
Foam:	Electrical:	Vinyl:
California Technical Bulletin 117 - Section E	100-240VAC 50/60Hz 5A	California T.B. 117-2013 NFPA 260- Cover Fabric - Class I UFAC Fabric - Class I

CERTIFICATIONS

Industry Standards:	Warranty:	Standard Color:
UL ANSI/AAMI ES60601-1 (2005) + AMD (2012), CAN/CSA C22.2 No. 60601 (2014)	2 Years	

PRODUCT BENEFITS

- Optional Trendelenburg to stabilize patient blood pressure
- Optional Fowler positioning to increase patient comfort
- Help protect sonographers from potential career-ending injuries by improving ergonomics
- Bariatric capable

OPTIONAL FEATURES

- Trendelenburg/Reverse Trendelenburg:
- 15°/15° Fowler (Head Up) Positioning: 65°
- Single Pedal Braking
- Safety Rails
- Sonographer Extension
- Foot Controller
- Safe-T-Wedge
- Pediatric/Geriatric Adapter
- IV Pole and Holder
- Padded Arm Board
- Memory Positioning