


GSPM Table™

The ultimate platform for General Surgery and Pain Management

MEDICAL POSITIONING, INC.
www.medicalpositioning.com
800-593-3246



STANDARD FEATURES

Height Adjustability:	Maximum Patient Capacity:	Imaging Area:	Trendelenburg/Reverse Trendelenburg:	Lateral Tilt:
31" to 41"	500 lbs	53"	20°/20°	20° Left and Right
Control System:	Braking:	Additional Features:		
Hand Controller	Individual Locking Casters	Cervical Access Shaped Top Prone Procedure Facial Opening	Carbon Fiber Top and 2" Pad Patient Safety Straps	C-arm Accessible Base

PRODUCT SPECIFICATIONS

Length:	Width:	Weight:
84"	22"	556 lbs
Foam:	Electrical:	
California Technical Bulletin 117 - Section E	120 VAC, 60 Hz, 1.6A Max	

PRODUCT BENEFITS

- Low attenuation reduces scatter for better imaging
- Cantilever design for better C-Arm accessibility during pain management procedures
- Bariatric capable

CERTIFICATIONS

Industry Standards:	Warranty:	Standard Color:
UL 60601-1, CAN/CSA C22.2 No. 601.1	5 Years	

OPTIONAL FEATURES

- Foot Controller
- Carbon Fiber Arm Boards