

# HUT Table™

The ultimate platform for syncope testing

MEDICAL POSITIONING, INC.  
www.medicalpositioning.com  
800-593-3246

6.5°/second tilt rate

2-way sonographer drop section/Patient back support (optional)

15° Trendelenburg and  
85° Reverse  
Trendelenburg

Weight capacity up to 350 lbs

Imaging drop section (optional)

Patient safety straps

31" to 39" height

Built-in folding foot platform

Hand controller  
with memory  
positioning



Optional features may be shown. Actual color may vary depending on monitor resolution and print quality.

## STANDARD FEATURES

<b>Height Adjustability:</b>	<b>Maximum Patient Capacity:</b>	<b>Tilt Rate:</b>	<b>Trendelenburg/Reverse Trendelenburg:</b>	<b>Control System:</b>
31" to 39"	350 lbs	6.5°/second	15°/85°	Hand Controller with Memory Positioning

<b>Braking:</b>	<b>Additional Features:</b>
Individual Locking Casters	Manual Lowering System Built-In Folding Foot Platform Patient Safety Straps

## PRODUCT BENEFITS

- Manual down feature for emergency situations
- Optional multi-modality features allow more product flexibility
- Bariatric capable folding foot platform
- Tilt rate meets ACC guidelines for syncope testing

## PRODUCT SPECIFICATIONS

<b>Length:</b>	<b>Width:</b>	<b>Weight:</b>
73"	30"	380 lbs
<b>Foam:</b>	<b>Electrical:</b>	<b>Vinyl:</b>
California Technical Bulletin 117 - Section E	120 VAC, 60 Hz, 1.6A Max 230 VAC, 50 Hz, 0.8A Max	California T.B. 117-2013 NFPA 260- Cover Fabric - Class I UFAC Fabric - Class I

## OPTIONAL FEATURES

- Drop Sections
  - Imaging Drop Section
  - 2-Way Sonographer Drop Section/Patient Back Support
- Sonographer Extension
- Safe-T-Wedge
- Pediatric/Geriatric
- Adapter IV Pole and Holder
- Head Extension
- Arm Board
- Safety Rails
- Foot Control

## CERTIFICATIONS

<b>Industry Standards:</b>	<b>Warranty:</b>	<b>Standard Color:</b>
UL 60601-1, CAN/CSA C22.2 No. 601.1	1 Year	