

RR HUT Table™

The ultimate platform for syncope and bradycardia testing

MEDICAL POSITIONING, INC.
www.medicalpositioning.com
800-593-3246




STANDARD FEATURES

Height Adjustability:	Maximum Patient Capacity:	Emergency Tilt Rate:	Trendelenburg/Reverse Trendelenburg:
31" to 39"	600 lbs	70° reverse Trendelenburg -level in 7 seconds (SWL for 600 lb patient)	15°/90°
Braking:	Additional Features:		PRODUCT BENEFITS
Individual Locking Casters	1 button emergency return to level Programmable home memory position Programmable tilt angle: 5° increments Built-in folding foot platform		
Control System:	Patient safety straps Battery backup		
Hand Controller with Memory Positioning and Battery Backup			7 second emergency response time when Rapid Response engaged Optional multi-modality features allows more product flexibility Bariatric capable folding foot platform Tilt rate meets ACC guidelines for syncope testing

PRODUCT SPECIFICATIONS

Length:	Width:	Weight:
73"	30"	380 lbs
Foam:	Electrical:	Vinyl:
California Technical Bulletin 117 - Section E	120 VAC, 60Hz, 6.5A 230 VAC, 50 Hz, 3.4A	California T.B. 117-2013 NFPA 260- Cover Fabric - Class I UFAC Fabric - Class I

CERTIFICATIONS

Industry Standards:	Warranty:	Standard Color:
UL 60601-1, CAN/CSA C22.2 No. 601.1	3 Years	

Latch
GUA-1304

OPTIONAL FEATURES

- Drop Sections
 - Imaging Drop Section
 - 2-Way Sonographer Drop Section/Patient Back Support
- Sonographer Extension
- Safe-T-Wedge
- Pediatric/Geriatric
- Adapter IV Pole and Holder
- Head Extension
- Arm Board
- Safety Handrails
- Foot Control