

RR HUT Table™

The ultimate platform for syncope and bradycardia testing

MEDICAL POSITIONING, INC.
www.medicalpositioning.com
800-593-3246



2-way sonographer drop section/Patient back support (optional)

Weight capacity up to 600 lbs

Patient safety straps

Built-in folding foot platform

Emergency tilt rate ability to bring table from 70°-level in 7 seconds (SWL for 600 lb patient)

Programmable tilt angle: 5° increments

Imaging drop section (optional)

15° Trendelenburg and 90° Reverse Trendelenburg

Battery backup

31" to 39" height

Hand controller with one button emergency return to level

Programmable home memory position


STANDARD FEATURES

Height Adjustability:	Maximum Patient Capacity:	Emergency Tilt Rate:	Trendelenburg/Reverse Trendelenburg:
31" to 39"	600 lbs	70° reverse Trendelenburg -level in 7 seconds (SWL for 600 lb patient)	15°/90°
Braking:	Additional Features:	PRODUCT BENEFITS	
Individual Locking Casters	1 button emergency return to level Programmable home memory position Programmable tilt angle: 5° increments Built-in folding foot platform Patient safety straps Battery backup	7 second emergency response time when Rapid Response engaged Optional multi-modality features allows more product flexibility Bariatric capable folding foot platform Tilt rate meets ACC guidelines for syncope testing	
Control System:			
Hand Controller with Memory Positioning and Battery Backup			

PRODUCT SPECIFICATIONS

Length:	Width:	Weight:
73"	30"	350 lbs
Foam:	Electrical:	Vinyl:
California Technical Bulletin 117 - Section E	120 VAC, 60 Hz, 6.0A Max	California T.B. 117-2013 NFPA 260- Cover Fabric - Class I UFAC Fabric - Class I

CERTIFICATIONS

Industry Standards:	Warranty:	Standard Color:
UL 60601-1, CAN/CSA C22.2 No. 601.1	3 Years	

OPTIONAL FEATURES

- Drop Sections
 - Imaging Drop Section
 - 2-Way Sonographer Drop Section/Patient Back Support
- Sonographer Extension
- Safe-T-Wedge
- Pediatric/Geriatric
- Adapter IV Pole and Holder
- Head Extension
- Arm Board
- Safety Handrails
- Foot Control