



Convertible platform for 360° access to breast

Chair to table positioning for multi use

The UltraDBI maximizes department space by incorporating upright and decubitus stereotactic breast biopsy, as well as mammography into one bariatric capable platform, increasing throughput and increasing continuity of care.


- Easy transition from chair to table provides an all-encompassing platform for varying patient population
- Accommodates both upright and decubitus positions for stereotactic breast biopsy, as well as mammography and breast ultrasound
- Bariatric capable patient weight capacity of 500 lbs establishes continuity of care
- Reduce turnover time and increase throughput, even in small exam rooms, with interchangeable head and foot rests eliminating the need to reposition the table to accommodate each patient
- Accessories Include: removable arm rests, back support, interchangeable, adjustable and removable head supports

STANDARD FEATURES

Height Adjustability:	Maximum Patient Capacity:	Trendelenburg:	Fowler:	Drop Sections:	Supported Patient Positions:
25" to 34.5"	500 lbs	15° to -15°	90°	1 on each side: one for physician and one for receptor	Upright in a chair and Decubitus on a table
Foot Platform:	Control System:	Braking:	Additional Features:		
Adjustable foot platform in the chair position	Hand Controller	Single pedal braking	Adjustable Arm Rest	Detachable Head Support	Adjustable Back Support Adjustable Head Support

PRODUCT SPECIFICATIONS

Length:	Surface Width:	Length of Wheel Base:	Foam:	Electrical:	Vinyl:
76"	28"	30"	California Technical Bulletin 117 - Section E	100-230 VAC, 50/60 Hz, Max 5 A	California T.B. 117-2013 NFPA 260- Cover Fabric - Class I UFAC Fabric - Class I

Certifications:	Warranty:	Standard Color:	OPTIONAL FEATURES	
UL ANSI/AAMI ES60601-1 (2005) + AMD (2012), CAN/CSA C22.2 No. 60601 (2014)	1 Year		IV Pole and Holder	Foot Control
			Battery	Underbed Lights
			Back Support Wedge Cushion	Instrument Tray