

Carotid head support (optional)

30° reverse Trendelenburg
and 15° Trendelenburg

Sonographer drop section/patient back
support and rapid release imaging drop
section (optional)

65° Fowler positioning

Bariatric capable patient
weight capacity up to 600 lbs.

Static load capacity
up to 1,000 lbs.

Lightweight/removable/
adjustable foot platform

Battery (optional)

Single pedal
braking (optional)

ADA compliant
19" to 38" at top
of mattress



Multi-modality platform



Transfer height of 19"
at top of mattress

ADA compliant 19" low height to 38" maximum height

The VasScan X helps minimize the risk of **repetitive stress injury**, **reduce exam time and improve results** while increasing **bariatric capability** and meeting the **ADA** height transfer requirement.

- 19"-38" height range meets the Americans with Disabilities Act ADA low height requirement of 17"-19" for patient transfer, while ensuring ergonomic standing posture for staff
- 30° reverse Trendelenburg increases venous pooling, prompting rapid recognition and accurate results
- The weight-bearing adjustable and removable foot platform accommodates patients of varying heights while ensuring patient safety during lower venous studies
- Bariatric capable patient weight capacity of 600 lbs establishes continuity of care
- 65° Fowler positioning increases comfort of/for patients with shortness of breath
- Optional imaging and sonographer drop sections provide unhindered access and reduce image acquisition time for echocardiograms

STANDARD FEATURES

Height Adjustability:	Maximum Patient Capacity:	Fowler (Head Up):	Trendelenburg/Reverse Trendelenburg:	Foot Platform:
19" to 38"	600 lbs Maximum Load Capacity: 1,000 lbs	65°	15°/30°	Manually adjustable and removable

Control System:	Braking:	Additional Features:
Hand Controller w/ Memory positioning	Individual Locking Casters	Hand Control Storage Foley Bag/Accessory Hook

PRODUCT SPECIFICATIONS

Length:	Width:
80.5"	28"

Foam:	Electrical:	Vinyl:
California Technical Bulletin 117 - Section E	100-240VAC 50/60Hz 5A	California T.B. 117-2013 NFPA 260- Cover Fabric - Class I UFAC Fabric - Class I

CERTIFICATIONS

Industry Standards:	Warranty:	Standard Color:
UL ANSI/AAMI ES60601-1 (2005) + AMD (2012), CAN/CSA C22.2 No. 60601 (2014)	2 Years	

NEW AND IMPROVED FEATURES

ADA Compliant transfer height of 19"	Underbed lights
Maximum height of 38" ergonomically accommodates standing staff	Battery operation - avoid cords when preferred
Increased Patient weight capacity to 600 lbs	Barrier free edges with new hideaway brackets
Increased warranty to 2 years	Updated and durable hand control
Memory Positioning - Store up to 2 positions with the touch of a button	Ability to utilize both hand and foot controls
	Carotid Head Support to ergonomically accommodate all scanning needs

OPTIONAL FEATURES

Echocardiogram Drop Sections	Safety handrails w/ hideaway brackets
-Sonographer Extension	Fold-down safety handrails
-Imaging Drop Section	Safe-T-Wedge
-2-Way Sonographer Drop Section/Patient Back Support	Pediatric/Geriatric Adapter
-Non-Pinch Flaps	IV Pole and Holder
Single Pedal Braking	Padded Arm Board
Foot Controller (non-memory)	Paper Roll Holder and Cutter
Battery Operation - up to a full day of operation with no need for plugging in	Carotid Head Support
	Underbed lights