

— UltraMamm™

The ultimate platform for Mammography and Breast Biopsy

MEDICAL POSITIONING, INC.
www.medicalpositioning.com
800-593-3246

Weight capacity up to 500 lbs.

90° Fowler Positioning

Articulating patient arm boards
and patient back supports

Battery operated

22" to 38" height

Electric built-in calf section
and folding foot platform




STANDARD FEATURES

Height Adjustability: 22" to 38"	Maximum Patient Capacity: 500 lbs	Fowler (Head Up): 0 - 90°	Calf Section: 0 - 83°	Braking: Individual Locking Casters
Control System: Hand Controller Included Battery Functionality	Foot Down Section: Electrically Adjustable Built-In Folding Foot Platform	Arm Boards: Dual Articulating Patient Arm Boards and Back Support		

PRODUCT SPECIFICATIONS

Length: 72.5"	Width: 28.5"	Weight: 350 lbs
Foam: TB 117-2013 Section E	Electrical: 100 - 240 VAC, 50/60 Hz, 4.5A Max	Vinyl: TB 117-2013 NFPA 260 Cover Fabric - Class 1 UFAC Fabric - Class 1

CERTIFICATIONS

Industry Standards: UL ANSI/AAMI ES60601-1 (2005) + AMD (2012), CAN/CSA C22.2 No. 60601 (2014)	Warranty: 2 Years	Standard Color: 
--	-----------------------------	---

PRODUCT BENEFITS

- Access lesions in upper/outer quadrant
- Compatible with any stereotactic add-on system
- Enhance ergonomic posture as technologists work upright
- Bariatric capable

OPTIONAL FEATURES

- Carotid Head Rest
- Foot Control
- Handrails with Hideaway Brackets
- Head Rail
- IV Pole and Holder
- Paper Roll Holder and Cutter
- Push Handles
- Safety Straps
- Single Pedal Braking
- Underbed LED Lighting