

MultiView™

BARIUM SWALLOW CHAIR

medicalpositioning.com
800-593-3246



FOLD-AWAY SEAT BACK NO
NEED TO REMOVE AND STORE

90° LEFT SEAT ROTATION

0° TO 75° CALF SECTION

90° FOWLER POSITIONING

600 LBS. WEIGHT CAPACITY



19" TRANSFER HEIGHT



LARGE PATIENT
TRANSPORT HANDLE

INTEGRATED FOLD-AWAY
SAFETY HANDRAILS

OXYGEN TANK HOLDER

HAND CONTROL WITH
PROGRAMMABLE OPTION

PATIENT SAFETY STRAPS

LARGE FOLD-AWAY
FOOT REST

BATTERY OPERATED

ELECTRIC LOCKING CASTERS

RADIOLOGY SERIES



STANDARD FEATURES*

Height Adjustability:	Load Capacity:	Fowler:
19"-30" ADA Compliant	600 lbs.	0°-90°
Calf Section:	Handrails:	Power Option:
0°-75°	Fold-Away	Battery and A/C
AP Radiolucent Imaging Window:	Seat Rotation:	Braking:
32"x13"	90° for AP positioning	Electric Locking

ADDITIONAL FEATURES

- | | |
|---------------------------|---------------------------------------|
| Accessory Hook | Patient Safety Straps |
| Electric Locking Casters | Integrated Fold-Away Safety Handrails |
| Fold-Away Seat Back | Large Patient Transport Handle |
| Large Fold-Away Foot Rest | Oxygen Tank Holder |

OPTIONAL FEATURES

- | | |
|----------------------------|------------------------|
| IV Pole and Holder | Retractable Power Cord |
| Charger and Second Battery | Head Safety Strap |

ADDITIONAL PRODUCT INFORMATION*

Table Length:	Control System:
72.3"	Hand Controller
Product Weight:	Standard Color:
335 lbs.	Grey
Warranty	
2 Years	

* All specifications and claims are subject to final testing approval

DIMENSIONS

