

MultiView™

STANDARD FEATURES

MedicalPositioning
800.593.3246
info@medicalpositioning.com

RADIOLOGY SERIES

MultiView™

BARIUM SWALLOW CHAIR

Height Adjustability: 19"-30" ADA Compliant	Load Capacity: 600 lbs.	Fowler: 0°-90°
Calf Section: 0°-75°	Handrails: Fold-Away	Power Option: Battery or A/C
AP Radiolucent Imaging Window: 32"x13"	Seat Rotation: 90° for AP positioning	Braking: Electric Locking Caster

ADDITIONAL FEATURES



Hand Controller



Accessory Hook



Electric Locking Casters



Fold-Away Seat Back



Patient Safety Straps



Integrated Fold-Away Safety Handrails



Large Patient Transport Handle




Oxygen Tank Holder



Large Fold-Away Foot Rest


OPTIONAL FEATURES



Charger and Second Battery



IV Pole and Holder



Tumble Forms 2 Deluxe Floor Sitter®
Small, Medium, Large

ADDITIONAL PRODUCT INFORMATION

Table Length: 72.3"	Weight: 335 lbs.
Width: 26"	Control System: Hand Controller
Warranty: 2 Years	Standard Color: Grey



INNOVATIVE AND VERSATILE

The MultiView™ is the future of modified barium swallow study platforms. It meets the most current ADA standard for height adjustability down to 19". This modified barium swallow study chair features fold-Away safety handrails, a fold-away self-storing back section, and a bariatric patient weight capacity of 600 lbs.* Our video fluoroscopy swallow study chair is an innovative, versatile, safe, and comfortable product.


MedicalPositioning

MEDICALPOSITIONING.COM

MPI maintains a policy of continual product improvement and reserves the right to change features, specifications, and prices without prior notification. Check with MPI for latest information. FORM MULTIVIEWBRO 110823 REV A Copyright © 2023 Medical Positioning, Inc. All rights reserved. Printed in the USA.

MultiView™ FEATURES

The MultiView™ promotes greater patient access with unmatched benefits including ADA compliant low transfer heights and bariatric patient support. With versatile options video swallow studies can be performed seamlessly while establishing continuity of care.



REDUCE COSTS AND INCREASE PROFITS

Reduce Patient Risk →
 1,000,000: Hospital falls per year extending stays by 6.3 days
 \$64,526: Average direct cost of a fall with injury

Reduce Staff Injury →
 25%: Worker's compensation claims caused by patient handling injuries
 \$15,600: Cost of each claim

Increase ADA Compliance →
 2,700,000: Americans rely on wheelchairs
 19": Wheelchair accessible. Easily increase ADA compliance, with a multi-functional chair

Increase Patient Satisfaction →
 High patient scores lead to bonuses for hospitals

Increase Staff Retention →
 20%: of RN's leave because of risk associated with their work
 \$52,350: to replace an RN

INCREASE PROFITS
 2008: Medicare does not cover the costs associated with inpatient falls. Avoid costs for workplace compensation claims and replacing staff