

MultiView™

BARIUM SWALLOW CHAIR

medicalpositioning.com
800-593-3246



FOLD-AWAY SEAT BACK NO
NEED TO REMOVE AND STORE

90° LEFT SEAT ROTATION

0° TO 75° CALF SECTION

90° FOWLER POSITIONING

600 LBS. WEIGHT CAPACITY



19" TRANSFER HEIGHT



LARGE PATIENT
TRANSPORT HANDLE

INTEGRATED FOLD-AWAY
SAFETY HANDRAILS

OXYGEN TANK HOLDER

HAND CONTROL WITH
PROGRAMMABLE OPTION

PATIENT SAFETY STRAPS

LARGE FOLD-AWAY
FOOT REST

BATTERY OPERATED

ELECTRIC LOCKING CASTERS

RADIOLOGY SERIES



STANDARD FEATURES

| | | |
|---------------------------------------|------------------------|----------------------|
| Height Adjustability: | Load Capacity: | Fowler: |
| 19"-30" ADA Compliant | 600 lbs. | 0°-90° |
| Calf Section: | Handrails: | Power Option: |
| 0°-75° | Fold-Away | Battery and A/C |
| AP Radiolucent Imaging Window: | Seat Rotation: | Braking: |
| 32"x13" | 90° for AP positioning | Electric Locking |

ADDITIONAL FEATURES

- | | |
|---------------------------|---------------------------------------|
| Accessory Hook | Patient Safety Straps |
| Electric Locking Casters | Integrated Fold-Away Safety Handrails |
| Fold-Away Seat Back | Large Patient Transport Handle |
| Large Fold-Away Foot Rest | Oxygen Tank Holder |

OPTIONAL FEATURES

- | | |
|-------------------------------------|------------------------|
| IV Pole and Holder | Retractable Power Cord |
| Charger and Second Battery | Head Safety Strap |
| Tumble Forms 2 Deluxe Floor Sitter® | |

ADDITIONAL PRODUCT INFORMATION

| | |
|------------------------|------------------------|
| Table Length: | Control System: |
| 72.3" | Hand Controller |
| Product Weight: | Standard Color: |
| 335 lbs. | Grey |
| Warranty | |
| 2 Years | |

DIMENSIONS

