

MultiView®

VIDEO SWALLOW STUDIES
MODEL: 5562

KEY STANDARD FEATURES:

- ✓ Fold-Away Seat Back Section, no need to remove and store
- ✓ 90° Patient's Left Seat Rotation
 - ✓ 32" x 13" AP Radiolucent Imaging Window
- ✓ Battery and A/C Power Option
- ✓ 19"-30" Height ADA Compliant
- ✓ 90° Fowler
- ✓ 75° Calf Section
- ✓ Integrated Folding Patient Foot Board
- ✓ Patient Safety Straps
- ✓ Large Patient Transport Handle
- ✓ Electric Casters
- ✓ 90° Fowler
- ✓ 600 lbs. Patient Weight Capacity
- ✓ 2 Year Warranty



ADA COMPLIANT

GSPM Table™

C-ARM IMAGING
MODEL: 60362

ADDITIONAL FEATURES:

- ✓ Carbon Fiber Dual Articulating Arm Boards
- ✓ Foot Control

KEY STANDARD FEATURES:

- ✓ C-Arm Accessible Cantilever Style Base
- ✓ 53" Carbon Fiber Imaging Surface
- ✓ Cervical Access Shaped Top with Prone Procedural Facial Opening
- ✓ 500 lbs. Patient Weight Capacity
- ✓ 31"-41" Height
- ✓ 5 Year Warranty
- ✓ +/-20° Lateral Roll
- ✓ +/-20° Trendelenburg



HÅG Capisco®

ADDITIONAL FEATURES:

- ✓ SHORT 16.25"-21.5" Height Column
- ✓ TALL 22.5"-32" Height Column
- ✓ Foot Base

KEY STANDARD FEATURES:

- ✓ 19"-26" Height
- ✓ Adjustable Seat Depth, Back Height and Tilt Resistance
- ✓ Four Direction Seating Support



ABOUT US

MedicalPositioning, founded in 1988 as American Echo, introduced the EchoBed® with ergonomic patient access windows developed to address clinician workflow and repetitive-motion considerations during sonography procedures. Today, MedicalPositioning offers ergonomic diagnostic imaging platforms that support clinician access and patient positioning across multiple specialties, including cardiology, vascular, radiology, breast biopsy, surgery, pain management, and speech pathology. Products are manufactured in Lenexa, Kansas.



800.593.3246 - INFO@MEDICALPOSITIONING.COM - 9732 PFLUMM RD - LENEXA, KS 66215



MedicalPositioning maintains a policy of continual product improvement and reserves the right to change features, specifications, and prices without prior notification. Check with MedicalPositioning for latest information. FORM MPIMINI 051526 REV A Copyright © 2026 MedicalPositioning, Inc. All rights reserved. Printed in the USA.

Proven



MedicalPositioning PRODUCT LINES



Innovative



Ergonomic



medicalpositioning.com

ADA COMPLIANT

EchoBed®

ADDITIONAL FEATURES:

85° Fowler

+/-15° Trendelenburg

Electric Patient Back Support Positioning

KEY STANDARD FEATURES:

Dual Rapid Release Imaging Drop Sections:

Series 700: 17" – 31" US Access Board Compliant

Series 900: 19"-37" ADA Compliant

Fold-Away Safety Handrails

Electric Locking Casters

750 lbs. Patient Weight Capacity

5 Year Warranty

CARDIAC IMAGING

MODEL: EB724, EB723, EB924, EB923

THE ALL NEW



SIMILAR MODEL

EchoTable®

MODEL: 251, 272, 273, 283



Stress EchoBed®

KEY STANDARD FEATURES:

Adjustable/Removable Ergometer

10 pre-programmed, 4 programmable, 1 manual protocol

Patient Coach with audible and visual patient guidance

Computer Controller/Interface

Dual Imaging Drop Sections:

Rapid Release Imaging Drop Section

2-Way Sonographer Drop Section/Patient Back Support

24"-34" Height

+/-15° Trendelenburg

40° Patient Left Lateral Roll

Four Point Safety System

500 lbs. Patient Weight Capacity

5 Year Warranty

CARDIAC IMAGING

MODEL: 1736



RR HUT®Table

ADDITIONAL FEATURES:

Dual Imaging Drop Sections:

Rapid Release Imaging Drop Section

2-Way Sonographer Drop Section/Patient Back Support

KEY STANDARD FEATURES:

One button emergency tilt rate 70° to level in 7 seconds (SWL for 600 lbs. patient)

31"-39" Height

15° Trendelenburg/90° Reverse Trendelenburg

Patient Safety Straps

Integrated Folding Patient Foot Board

Battery Backup

600 lbs. Patient Weight Capacity

3 Year Warranty

CARDIAC IMAGING

TILT STUDIES

VASCULAR IMAGING

MODEL: 1056, 1256



SIMILAR MODEL

HUT®Table

MODEL: 1031, 1231



Cath Ergometer™

ADDITIONAL FEATURES:

Storage Cart

KEY STANDARD FEATURES:

Custom Fit Adjustable/Removable Ergometer to Specific Cath or Nuclear Table of choice

10 pre-programmed, 4 programmable, 1 manual protocol

Patient Coach with audible and visual patient guidance

Computer Controller/Interface

2 Year Warranty

CARDIAC IMAGING



UltraMamm®

ADDITIONAL FEATURES:

Integrated Dual Articulating Arm Boards

+/-15° Trendelenburg

Single Pedal Braking

Removable Safety Handrails w/Hide-Away Brackets

Patient Safety Strap

KEY STANDARD FEATURES:

22"-38" Height

90° Fowler

83° Calf Section

Integrated Folding Patient Foot Board

Battery and A/C Power Options

500 lbs. Patient Weight Capacity

2 Year Warranty

DIAGNOSTIC IMAGING

BREAST IMAGING

MODEL: 052, 052-T, 952, 952-T



SIMILAR MODEL

Ultra DBI™

MODEL: 9400



VasScan® Table X

ADDITIONAL FEATURES:

Dual Imaging Drop Sections:

Rapid Release Imaging Drop Section

2-Way Sonographer Drop Section/Patient Back Support

Single Pedal Braking

Safety Handrails

KEY STANDARD FEATURES:

19"-38" Height ADA Compliant

65° Fowler

15° Trendelenburg/30° Reverse Trendelenburg

Adjustable/Removable Patient Foot Board

600 lbs. Patient Weight Capacity

2 Year Warranty

CARDIAC IMAGING

VASCULAR IMAGING

MODEL: 084, 284, 984



UltraScan Versa®

ADDITIONAL FEATURES:

Dual Imaging Drop Sections:

Rapid Release Imaging Drop Section OR

2-Way Sonographer Drop Section/Patient Back Support

+/-15° Trendelenburg or 15° Trendelenburg/25° Reverse Trendelenburg

Integrated Folding Patient Foot Board

Removable Safety Handrails w/Hide-Away Bracket OR Patient Arm Rests

Patient Safety Strap

Single Pedal Braking

KEY STANDARD FEATURES:

Adjustable/Hideaway Stirrups

22"-38" Height

90° Fowler

83° Calf Section

500 lbs. Patient Weight Capacity

2 Year Warranty

DIAGNOSTIC IMAGING

BREAST IMAGING

CARDIAC IMAGING

VASCULAR IMAGING

MODEL: 087, 287, 987, 095, 295, 995, 092



UltraScan Versa Premier®

KEY STANDARD FEATURES:

22"-38" Height

Dual Imaging Drop Sections:

Integrated Dual Articulating Arm Boards

OR Dual Imaging Drop Sections

Adjustable/Hideaway Stirrups

90° Fowler

83° Calf Section

15° Trendelenburg/25° Reverse Trendelenburg

Integrated Folding Patient Foot Board

Single Pedal Braking

Removable Safety Handrails w/Hide-Away Brackets

500 lbs. Patient Weight Capacity

5 Year Warranty

DIAGNOSTIC IMAGING

BREAST IMAGING

CARDIAC IMAGING

VASCULAR IMAGING

MODEL: 292, 992



MEDICALPOSITIONING.COM